


# BioNgraft<sup>TM</sup> Flex

Flexible Synthetic Bone Substitute



## Product Features

### Easy to Use

Biongraft Flex provides great handling with high elasticity for specific cases.

### 100% Synthetic

Contains no tissue of human or animal origin therefore carries no risk of disease transmission.

### Bioresorbable

With its optimized porous structure and chemical composition, Biongraft Flex is suitable for the continuous remodeling cycle of healthy bone. b-TCP resorbs over time and be replaced with bone during the healing process.

### Osteoconductive

Act as a scaffold and support bone tissue regeneration. Similar to the mineral found in bone tissue.

Highly Interconnected porosity with excellent mechanical resistance.

# BioNgraft<sup>TM</sup> Flex

Flexible Synthetic Bone Substitute

## Composition

Biongraft Flex is composed of b-TCP and resorbable PLA based Synthetic Polymer. The polymer of the Flex is metabolite lactic acid and is degraded to CO<sub>2</sub> and H<sub>2</sub>O.

## Indications

Biongraft Flex is indicated to be gently packed or placed into bony voids or gaps of the skeletal system (i.e. the extremities spine and pelvis) and to fill cavities or defects resulting from cysts, tumors, or other causes major disc collapse.

## Easy Application

- Biongraft Flex can be applied directly or in combination with bone marrow aspirate/blood to the surgical site.
- Wetting Biongraft Flex increases flexibility.
- Biongraft Flex can be cut to fit into a lumbar or cervical cage.

## Specifications

References	Description	Size	Packaging
<b>TBGF25254</b>	Flexible Bone Strip	25x25x4mm	1 Unit/Box
<b>TBGF701104</b>	Flexible Bone Strip	70x110x4mm	1 Unit/Box
<b>TBGF25804</b>	Flexible Bone Strip	25x80x4mm	1 Unit/Box
<b>TBGF25805</b>	Flexible Bone Strip	25x80x5mm	1 Unit/Box
<b>TBGF25806</b>	Flexible Bone Strip	25x80x6mm	1 Unit/Box
<b>TBGF60504</b>	Flexible Bone Strip	60x50x4mm	1 Unit/Box
<b>TBGF60505</b>	Flexible Bone Strip	60x50x5mm	1 Unit/Box
<b>TBGF60606</b>	Flexible Bone Strip	60x60x6mm	1 Unit/Box
<b>TBGF45806</b>	Flexible Bone Strip	45x80x6mm	1 Unit/Box
<b>TBGF35605</b>	Flexible Bone Strip	35x60x5mm	1 Unit/Box
<b>TBGF35606</b>	Flexible Bone Strip	35x60x6mm	1 Unit/Box
<b>TBGF75405</b>	Flexible Bone Strip	75x40x5mm	1 Unit/Box
<b>TBGF30306</b>	Flexible Bone Strip	30x30x6mm	1 Unit/Box

Produced exclusively for



10, Kassiopis Str., 17237,  
Ymittos, Greece  
Email: [info@eosmed.eu](mailto:info@eosmed.eu)  
Website: [www.eosmed.eu](http://www.eosmed.eu)  
Tel: +302160033602

